

Columbia Coatings

# Technical Data Sheet

Columbia Coatings  
1173 Industrial Park Rd  
Columbia, TN 38401  
931-388-7730  
931-388-5573 FAX



<b>Number:</b>	THERMO-STRIP
<b>Name:</b>	Thermo Strip 2048 Liquid
<b>Type:</b>	Liquid Stripper

## Chemical Properties

<b>Appearance</b>	Clear amber liquid
<b>Odor</b>	Mild
<b>pH</b>	11

## Recommended Application

<b>Step #1</b>	Prepare a strip tank with lid and proper ventilation.
<b>Step #2</b>	Immerse part to be stripped in Thermo Strip Liquid.
<b>Step #3</b>	Heat tank to 165-180°F. <i>Stir solution with mixer during stripping process. Stripping time is typically 20-30 minutes.</i>
<b>Step #4</b>	Rinse parts using fresh water and pressure washer.

## Benefits

- Works well on coated carbon steel, stainless steel and aluminum
- Minimal airborne health hazard
- Low vapor pressure reduces VOC's, Environmentally friendly

## Health and Safety

- Please refer to SDS for complete health and safety information.

## Storage

- Store in a dry location at room temperature. Keep container closed and maintain all original markings and labels. Keep this container and vapors from this container away from heat, sparks, flame, and other ignition sources. Recommended storage in ferrous, stainless steel and plastic containers only. Note: do not store in aluminum containers at any time.

[www.columbiacoatings.com](http://www.columbiacoatings.com)